Course Software Version 2014 November 2014 Edition Part Number 325511D-01

Copyright

© 2004-2014 National Instruments. All rights reserved.

LabVIEW[™] Core 3 Exercises

Under the copyright laws, this publication may not be reproduced or transmitted in any form, electronic or mechanical, including photocopying, recording, storing in an information retrieval system, or translating, in whole or in part, without the prior written consent of National Instruments Corporation.

National Instruments respects the intellectual property of others, and we ask our users to do the same. NI software is protected by copyright and other intellectual property laws. Where NI software may be used to reproduce software or other materials belonging to others, you may use NI software only to reproduce materials that you may reproduce in accordance with the terms of any applicable license or other legal restriction.

End-User License Agreements and Third-Party Legal Notices

You can find end-user license agreements (EULAs) and third-party legal notices in the following locations:

- + EULAs are located in the <National Instruments> $\SharedMDF\Legal\license directory.$
- $\label{eq:logal} equation and logal information in installers built with NI products.$

Trademarks

Refer to the NI Trademarks and Logo Guidelines at ni.com/trademarks for more information on National Instruments trademarks.

ARM, Keil, and µVision are trademarks or registered of ARM Ltd or its subsidiaries.

LEGO, the LEGO logo, WEDO, and MINDSTORMS are trademarks of the LEGO Group.

TETRIX by Pitsco is a trademark of Pitsco, Inc.

FIELDBUS FOUNDATION[™] and FOUNDATION[™] are trademarks of the Fieldbus Foundation.

EtherCAT[®] is a registered trademark of and licensed by Beckhoff Automation GmbH.

CANopen® is a registered Community Trademark of CAN in Automation e.V.

DeviceNet[™] and EtherNet/IP[™] are trademarks of ODVA.

Go!, SensorDAQ, and Vernier are registered trademarks of Vernier Software & Technology. Vernier Software & Technology and vernier.com are trademarks or trade dress.

Xilinx is the registered trademark of Xilinx, Inc.

Taptite and Trilobular are registered trademarks of Research Engineering & Manufacturing Inc.

FireWire[®] is the registered trademark of Apple Inc.

Linux® is the registered trademark of Linus Torvalds in the U.S. and other countries.

Handle Graphics[®], MATLAB[®], Real-Time Workshop[®], Simulink[®], Stateflow[®], and xPC TargetBox[®] are registered trademarks, and TargetBox[™] and Target Language Compiler[™] are trademarks of The MathWorks, Inc.

Tektronix[®], Tek, and Tektronix, Enabling Technology are registered trademarks of Tektronix, Inc.

The Bluetooth® word mark is a registered trademark owned by the Bluetooth SIG, Inc.

The ExpressCard[™] word mark and logos are owned by PCMCIA and any use of such marks by National Instruments is under license.

The mark LabWindows is used under a license from Microsoft Corporation. Windows is a registered trademark of Microsoft Corporation in the United States and other countries.

Other product and company names mentioned herein are trademarks or trade names of their respective companies.

Members of the National Instruments Alliance Partner Program are business entities independent from National Instruments and have no agency, partnership, or joint-venture relationship with National Instruments.

Patents

For patents covering National Instruments products/technology, refer to the appropriate location: Help»Patents in your software, the patents.txt file on your media, or the National Instruments Patent Notice at ni.com/patents.

Worldwide Technical Support and Product Information $\mathtt{ni.com}$

Worldwide Offices

Visit ni.com/niglobal to access the branch office websites, which provide up-to-date contact information, support phone numbers, email addresses, and current events.

National Instruments Corporate Headquarters

11500 North Mopac Expressway Austin, Texas 78759-3504 USA Tel: 512 683 0100

To comment on National Instruments documentation, refer to the National Instruments website at ni.com/info and enter the Info Code feedback.



Contents

Student Guide

	A. NI Certification	vii
]	B. Course Description	viii
(C. What You Need to Get Started	ix
]	D. Installing the Course Software	ix
]	E Course Goal	ix

Lesson 1

Developing Successful Applications Exercise 1-1 Review a Requirements Document. Exercise 1-2 Review User Stories.

Lesson 2

Organizing the Project			
Exercise 2-1	Create a Project Library	 	
Exercise 2-2	Resolve Project Conflicts		

Lesson 3

Creating an Applica	tion Architecture		
Exercise 3-1	Create a Queued Message Handler	 	 3-3
Exercise 3-2	Handshaking with Notifiers		
	8		

Lesson 4

Customizing the User Interface

Exercise 4-1	Create a User Interface Prototype	
Exercise 4-2	Create User Documentation	
Exercise 4-3	Initialize an Application from a File	
Exercise 4-4	Improve Application Usability	
	I II II J	

Contents

Lesson 5		
Managing and Loggi	ing Errors	
Exercise 5-1	Manage Errors	
Exercise 5-2	Log Errors	5-17
Lesson 6 Creating Modular Co	ode	
Exercise 6-1	Implement and Test a Data Logging Module	
Exercise 6-2	Integrate and Test Modules from the Top-Down	
Exercise 6-3	Integrate and Test Modules from the Bottom-Up	
Exercise 6-4	Implement Fuel Control Valve Functionality	
Exercise 6-5	Integrate and Sandwich Test a Module	
Annendix A		

Appendix A Additional Information and Resources

Student Guide

Thank you for purchasing the *LabVIEW Core 3* course kit. This course manual and the accompanying software are used in the three-day, hands-on *LabVIEW Core 3* course.


You can apply the full purchase of this course kit toward the corresponding course registration fee if you register within 90 days of purchasing the kit. Visit ni.com/training for online course schedules, syllabi, training centers, and class registration.

Note For course and exercise manual updates and corrections, refer to ni.com/info and enter the Info Code core3.

A. NI Certification

E

The *LabVIEW Core 3* course is part of a series of courses designed to build your proficiency with LabVIEW and help you prepare for NI LabVIEW certification exams. The following illustration shows the courses that are part of the LabVIEW training series. Refer to ni.com/training for more information about NI Certification.



B. Course Description

LabVIEW Core 3 introduces you to structured practices to design, implement, document, and test LabVIEW applications. This course focuses on developing hierarchical applications that are scalable, readable, and maintainable. The processes and techniques covered in this course help reduce development time and improve application stability. By incorporating these design practices early in your development, you avoid unnecessary application redesign, increase VI reuse, and minimize maintenance costs.

This course assumes that you have taken the LabVIEW Core 1 and LabVIEW Core 2 courses or have equivalent experience.

This course kit is designed to be completed in sequence. The course and exercise manuals are divided into lessons, described as follows.

In the course manual, each lesson consists of the following;

- An introduction that describes the purpose of the lesson and what you will learn
- A discussion of the topics in the lesson
- A summary quiz that tests and reinforces important concepts and skills taught in the lesson

In the exercise manual, each lesson consists of the following:

- A set of exercises to reinforce the topics in the lesson
- · Some lessons include optional and challenge exercise sections or additional exercises to complete if time permits

Note The exercises in this course are cumulative and lead toward developing a final application at the end of the course. If you skip an exercise, use the solution VI for that exercise, available in the <Solutions>\LabVIEW Core 3 directory, in later exercises.

C. What You Need to Get Started

Before you use this course manual, make sure you have the following items:

- □ Windows XP or later installed on your computer
- □ LabVIEW Professional Development System 2012 or later
- □ *LabVIEW Core* 3 course CD, containing the following folders:

Directory			Description
Exercises	F	Folder containing VIs	and other files used in the course
Solutions	F	Folder containing con	pleted course exercises

D. Installing the Course Software

Complete the following steps to install the course software.

- 1. Insert the course CD in your computer.
- 2. Follow the prompts to install the course material.


The installer places the Exercises and Solutions folders at the top level of the root directory. Exercise files are located in the <Exercises>\LabVIEW Core 3 directory.

Tip Folder names in angle brackets, such as <Exercises>, refer to folders in the root directory of your computer.

E. Course Goal

Given a requirements document for a LabVIEW development project, you will follow a software development process to design, implement, document and test the key application features in a manner that satisfies requirements for readability, scalability, and maintainability.







Exercise 5-1 Manage Errors

Goal

Modify your code to gracefully handle both critical and non-critical errors.

Scenario

In the sprint planning meeting for this iteration, the product owner chose to implement the following user stories:

- As a boiler operator, I want the boiler controller to prevent the boiler from entering an unsafe state as a result of a critical error.
- As a boiler operator, I want the system to be robust enough that minor errors do not result in system shut down, so that the boiler only shuts down if the error prevents safe operation.

Critical errors require shutting down the application because they cannot be otherwise managed. If they are not handled, critical errors may result in injury to personnel or damage to equipment. They may also place the system into an unmanageable state.

Non-critical errors are more akin to inconveniences that should not prevent continued execution of the application.

After discussing these user stories with the team, you identified and classified the following errors:

Table 5-1. Errors

Cause of Error	Critical or Non-Critical
Cancelling the Open File dialog when loading the INI file	Non-critical
The user chooses an invalid INI file	Critical
During development, you accidentally mistype a message string or case name	Critical

To address these errors, your team identifies the following tasks:


1. Modify Read Configuration Data.vi to handle cancellation of the Open File dialog.

- 2. Modify Read UI Constants.vi to handle selection of an invalid INI file.
- 3. Modify the Default case of the Message Handling Loop and Boiler Controller to handle typos when sending messages.
- 4. Modify the Message Handling Loop and Boiler Controller Loop to handle shutting down in the event of a critical error.

Design

The team decides on the following strategies for approaching each task:


- 1. Modify Read Configuration Data.vi to handle cancellation of the Open File dialog.
 - Modify Read Configuration Data.vi to load a default INI file if the user cancels the Open File dialog.
- 2. Modify Read UI Constants.vi to handle selection of an invalid INI file.
 - Modify Read UI Constants.vi to check that each key was properly read from the INI file.
 - Use an error ring to generate a custom error to indicate that the INI file is not valid.
- 3. Modify the Default case of the Message Handling Loop and Boiler Controller to handle typos when sending messages.
 - Use an error ring to generate a custom error to indicate that an invalid message was received by that loop.
- 4. Modify the Message Handling Loop and Boiler Controller Loop to handle shutting down in the event of a critical error.
 - Create a separate Error case in the Message Handling Loop and Boiler Controller Loop.
 - Modify Dequeue Message.vi to send you to the Error case instead of the Exit case in the event of an error.
 - Critical errors should be passed as message data from the bottom up (Boiler»Boiler Controller»Message Handling Loop)
 - The Error case of the Message Handling Loop sends you to the Initiate Stop case, as if you clicked Emergency Stop.



Implementation

Non-Critical Errors

- 1. Modify Read Configuration Data.vi to handle cancellation of the Open File dialog.
 - Run Main.vi. Cancel the Open File dialog where you would normally select the INI file. Notice that the File Dialog express VI generates Error 43 to indicate that the operation was cancelled. This error causes the application to exit. This is *not* a desirable behavior.
 - □ Modify Read Configuration Data.vi to load a default INI file if the user cancels the Open File dialog as shown in Figure 5-2.






Figure 5-2. Read Configuration Data VI Block Diagram Load Default INI File

- 1
- 2 Path constant
- Wire the **cancelled** output of the File Dialog express VI to the case selector. 3
- 2. Test the VI.
 - Run Main.vi and cancel the Open File dialog. Notice that the code continues executing normally, using the default file.

Critical Errors

- 1. Modify Read UI Constants to handle selection of an invalid INI file.
 - Run Main.vi. Select any file other than the INI file that you created. Notice that the range of Fuel Control Valve updated to 0-1 instead of 10-75. This did not generate an error, but the application cannot function without proper configuration of the Fuel Control Valve limits.
 - □ Modify Read UI Constants.vi as shown in Figure 5-3 to check that each key was properly read from the INI file.




Figure 5-3. Read UI Constants Using Error Ring

1 Error Ring—Generate a custom error to indicate that the INI file is not valid. This error shuts down the system.

□ Configure the Error ring with custom error code 5000 and the message Could not read UI constants from INI file.


- 2. Modify the Default case of the Message Handling Loop and Boiler Controller to handle typos when sending messages.
 - □ Modify the Reset event case of the Event Handling Loop to send the incorrect message, Resets, to the Message Handling Loop.





Change the Message string constant input to Resets. 1

Open the Default case of the Message Handling Loop. Notice that the code is already generating an error. Replace the Format into String and Error Cluster From Error Code functions with an error ring, as shown in Figure 5-5. This method of generating an error allows greater flexibility in how you generate the error.



1 Error Ring—Configure the error ring with Custom Error Code 5003 and the message The invalid message string "%s" was received in the Message Handling Loop.

3. Modify the Default case of the Boiler Controller.vi as shown in Figure 5-6.




Figure 5-6. Boiler Controller VI Default Case with Error Ring

- 1 Error ring—Configure the error ring with custom error code 5004 and the message The invalid message string "%s" was received in the Boiler Controller Loop.
 - Run Main.vi. Click **Reset**. Notice that the application hangs because Dequeue Message.vi sends the Message Handling Loop directly into the Exit case which causes the Message Handling Loop and Event Handling Loop stop, but no message is sent to the boiler controller to stop it and the boiler.
 - \Box To halt the application, select the Boiler window and press <Ctrl -.>.

- 4. Modify the Message Handling Loop and Boiler Controller Loop to handle shutting down in the event of a critical error.
 - □ Create an Error case in the Message Handling Loop and Boiler Controller Loop.





Figure 5-7. Message Handling Loop Error Case



5. Modify Dequeue Message.vi, as shown in Figure 5-9 to send you to the Error case instead of the Exit case in the event of an error.

Note You must close all open VIs and modify Dequeue Message.vi from the LabVIEW project. This VI is reentrant, so you must ensure that there is only one copy in memory before you can edit it.




Figure 5-9. Dequeue Message Block Diagram

- 2 Update the Message to Stop the Handler to Error.
- 3 Error cluster constant

1

Critical errors should be passed as message data from the bottom up—Boiler»Boiler Controller»Message Handling Loop.



Note Boiler Controller Error sends a message to the Message Handling Loop Error. Errors in the boiler are passed up to the Boiler Controller in the same way.

The Error cases should *not* wire the error cluster into the Enqueue Message. For the Boiler Controller, wire an error cluster constant to the shift register for each loop to avoid an endless error cycle. (The desired effect for the error will happen by launching the shutdown process.)

□ The Error case of the Message Handling Loop sends you to the Initiate Stop case, as if you clicked Emergency Stop.

Test the Application

- 1. Test the VI.
 - $\hfill\square$ Run the VI.
 - \Box Verify that:
 - Cancelling the Open File dialog does not stop the application.
 - Specifying an invalid INI file stops the application.
 - All previously-implemented functionality except for Reset still works.
 - Clicking **Reset** generates a critical error that halts the application.
- 2. Test the build specification.
 - \Box Build the executable from the build specification.
 - $\hfill\square$ Run Boiler Controller.exe.
 - \Box Verify the behavior from Step 1.
 - Cancelling the Open File dialog fails because the relative path to the default INI file is different for the EXE—in the data folder, instead of at the same level.

3. Modify Read Configuration Data.vi, as shown in Figure 5-10, to use a different relative path for EXEs.




Figure 5-10. Modify Read Configuration Data Block Diagram

- 4. Fix the forced error.
 - □ Update the Reset case of the Event Handling Loop to call Reset instead of Resets.

End of Exercise 5-1

Exercise 5-2 Log Errors

Goal

Log information about any error that causes a system shutdown.

Scenario

In the sprint planning meeting for this iteration, the product owner chose to implement the following user story:

• As a boiler operator, I want information about critical errors to be logged to disk, so that I can identify and troubleshoot the cause of the error at a later date.

As it is currently implemented, critical errors cause the application to shut down without providing any additional information to the end-user.

To implement this user story, your team identifies the following tasks:

- 1. Create the error log file in the same directory as the main application.
- 2. When a critical error occurs, log the error code and source.
- 3. If you cannot write to the application directory, display the error information to the user.
- 4. After logging the error, close the error log file reference.

Design

The team decides on the following strategies for approaching each task:

- 1. Create the error log file in the same directory as the main application.
 - The error log should only be created if a critical error actually occurs.
 - Create the error file in the Error case of the Message Handling Loop.
- 2. When a critical error occurs, log the current time, error code and source.
 - Pass the error cluster as Message Data when you send the Error message.
 - Extract the error code and source from the error cluster.
 - Create a subVI that writes the time, error code, and source to the error log file.

- 3. If you cannot write to the application directory, display the error information to the user.
 - If the error log reference is invalid, then you cannot write to the file.
 - Instead, launch a one-button dialog and display the information to the user.
- 4. After logging the error, close the error log file reference.

Implementation

- 1. Add a VI to write the time, error code, and source to the error log file.
 - □ Create a virtual folder called Logging in Support VIs folder the Boiler Controller Project Explorer window.
 - □ Navigate to <Exercises>\LabVIEW Core 3\External\Support VIs\Logging and copy Log Error.vi into the <Exercises>\ LabVIEW Core 3\Course Project\support\Logging directory.
 - □ Add Log Error.vi to the Logging virtual folder of the Boiler Controller Project Explorer window.

- 2. If you cannot write to the application directory, display the error information to the user.
 - □ Modify the error file in the Error case of the Message Handling Loop as shown in Figure 5-11.




Figure 5-11. Message Handling Loop Error Case

1

2

- 3. When a critical error occurs, log the current time, error code and source.
 - □ Modify the Error case of Boiler Controller.vi, as shown in Figure 5-12. You use the custom errors 5001 and 5002 in future exercises.




Figure 5-12. Boiler Controller Error Case

Test the Application

- 1. Test the VI.
 - □ Run the VI. Verify the following:
 - Specifying an invalid INI file stops the application and creates Error Log.txt.
 - Error Log.txt contains time, error code, and error source information similar to Figure 5-13.

Figure 5-13. Error Log.txt

Error Log.txt - Notepad	D	
<u>File Edit Format View H</u> elp		
4:14:54 PM 9/2/2014 5000 Read UI Constants.vi <err></err>		*
Could not read UI constants from INI file.		
		-

- □ Modify the Reset event case of the Event Handling Loop to send a **Resets** message to the Message Handling Loop.
 - Run the VI again.
 - Clicking Reset halts the application with a critical error.
 - Open Error Log.txt and view the new entry to the file.
 - Close the file.
 - Fix the Reset case to send a **Reset** message.
- \Box Verify the behavior when the application cannot write to the log file.
 - Set Error Log.txt to read-only.
 - Run the VI.

- Specify an invalid INI file. A dialog appears indicating the time, error code, and source.
- Set Error Log.txt to be writable.
- 2. Test the build specification.
 - □ Build the executable from the build specification.
 - □ Run Boiler Controller.exe.
 - \Box Verify the behavior from Step 1.
 - Specify an invalid INI file. You do not need to change the Reset case.

End of Exercise 5-2