

NATIONAL INSTRUMENTS™

NI Academic Site License

Installation Instructions

Thank you for purchasing National Instruments software. Please read these documents to familiarize yourself with the contents of your license.

Installation Order

National Instruments recommends that you install your NI Academic Site License software, as necessary, in the following order:

1. Software from the LabVIEW Core DVDs.
2. Software from the Extended Development Suite DVD.
3. Software from other Academic Site License options: Signal Processing and Communications Option, Control & Embedded Systems Option, and Circuit Design and Simulation Software.
4. NI Device Drivers.

Serial Numbers

Use the serial numbers included in your package to activate your software, as follows:

Department, College or Campus Licenses

If you purchased one of the following Academic Site Licenses, you will receive a single serial number on the Certificate of Ownership card. Use this serial number when installing and activating all software in your package.

- Department Teaching License
- College Teaching License
- Campus Teaching License

Department Teaching License with Options

If you purchased the Department Teaching License with individual options, such as the Control and Embedded Option, each DVD will contain a Certificate of Ownership card with a unique serial number. Use each serial number to activate software from its respective DVD.

Product Installation and Activation

Keep in mind the following guidelines when installing your NI software.

A. Single computer

You can install the software on individual computers. In this case, you need to activate the software on each computer.

You can choose from several ways to activate your software. When you launch the product after installation, a dialog box alerts you that activation is required to use the release version of the software, as well as the steps necessary for activation. If you cannot activate immediately, the software runs in evaluation mode. For more information on software activation, refer to the activation instructions included with the product, or visit ni.com/activate.

B. Multiple computers

If you have a large number of computers you would like to install the software on, you can make and distribute an image of the software.

In order to license software on the imaged computers, you can choose from two recommended methods:

1. NI Volume License Manager (NI VLM)
2. A computer image that contains an activated license file

For a description of both methods, refer to the National Instruments Web site at ni.com/info and enter the info code `rdhc15`.

NI VLM makes managing multiple licenses easy. For more information on NI VLM, contact your NI representative or refer to the National Instruments Web site at ni.com/vlp. To download NI VLM for FREE, refer to the National Instruments Web site at ni.com/info and enter the info code `NIVLMDOWNLOAD`.

Product Registration

Fill out the included registration card and return it to National Instruments. NI uses the contact information that you provide to facilitate technical support and to ensure you receive important updates.

Documentation

Documentation, including online help, for LabVIEW, LabWindows™/CVI™, Measurement Studio, and add-on software is installed with each product. Also, you can download product manuals from ni.com/manuals at no cost.

National Instruments, NI, ni.com, and LabVIEW are trademarks of National Instruments Corporation. Refer to the *Terms of Use* section on ni.com/legal for more information about National Instruments trademarks. Other product and company names mentioned herein are trademarks or trade names of their respective companies. For patents covering National Instruments products/technology, refer to the appropriate location: **Help»Patents** in your software, the `patents.txt` file on your media, or the *National Instruments Patent Notice* at ni.com/patents.